

**Roll No.16010421063**

**Batch No.G3**

**Name:Arya Nair**

### **Experiment 4: COD determination of waste-water**

**Objective:**

To determine chemical oxygen demand (COD) of water samples

**Chemical Oxygen Demand (COD):**

**Theory:**

COD is used as a measure of oxygen equivalent to organic matter content of a sample that is susceptible to oxidation by a strong chemical oxidant. For samples from a specific source, COD can be related empirically to BOD. COD determination has advantage over BOD determination in that the result can be obtained in about 5 hours as compared to 5 days required for BOD test.

The organic matter gets oxidized completely by  $K_2Cr_2O_7$  in the presence of  $H_2SO_4$  to produce  $CO_2$  and  $H_2O$ . The excess of  $K_2Cr_2O_7$  remained after the reaction is titrated with ferrous ammonium sulphate. The dichromate consumed gives the  $O_2$  required for oxidation of organic matter.

**Procedure:**

**Determination of COD of water sample**

1. Select the water sample.
2. To reflux the contents in the RB flask click the "switch on mantle" button.
3. Click "start titration" to titrate the contents.
4. Select the normality of ferrous ammonium sulphate (FAS).
5. Start titration & note the volume of titrant consumed when colour changes from bluish green to wine red. (Let the volume of titrant be  $V_2$  mL).
6. Repeat the same with the blank (Let the volume of the titrant be  $V_1$  mL).
7. COD calculated using the equation.

$$COD \left( \frac{mg}{L} \right) = \frac{Vol.FAS \text{ for sample} \times Normality \text{ of FAS} \times 8000}{Vol. \text{ of sample}}$$

### Observations and Calculations:

	Sample Type	Vol. of sample (mL)	Burette Reading (mL)		Vol. of FAS (mL)
			Initial	Final	
<b>Sample (3 h)</b>	Well water	20	0.0	4.8	4.8
	Tap	20	0.0	4.9	4.9
	Domestic	20	0.0	2.5	2.5
<b>Blank (0 min)</b>	Well water	20	0.0	5	5
	Tap	20	0.0	5	5
	Domestic	20	0.0	5	5

Water Analysis : Chemical Content

VARIABLES

Select test:  
COD

Titrant :  
Blank (20ml)

Normality : 0.5

Speed of titrant : 0.2

START  
RESET

RESULT

Titrant : Ferrous ammonium sulphate  
Titrant used : 5 ml  
Indicator : Ferroin

### COD of well-water

Potassium dichromate consumed =  $(V_{\text{blank}} - V_{\text{sample}}) = 0.2$

Normality of FAS = .....0.5....N.


Volume of the water sample = .....20.....mL.

$$\frac{\text{Vol. FAS for sample} \times \text{Normality of FAS} \times 8000}{\text{Vol. of sample}}$$

Therefore COD of the water sample =

= ...40.....ppm

Water Analysis : Chemical Content



Copyright © Amrita University 2009 - 2015

VARIABLES

**Select test:**  
COD

**Titrant:**  
Well water (20ml)

**Normality:** 0.5

**Speed of titrant:** 0.2

START  
RESET

RESULT

**Titrant:** Ferrous ammonium sulphate

**Titrant used:** 4.8 ml

**Indicator:** Ferroin

## COD of Domestic-water

Volume of FAS used=  $(V_2 - V_1)$  =.....2.5.....mL.

Normality of FAS =.....0.5.....N.


Volume of the water sample =.....20.....mL.

$$\frac{\text{Vol.FAS for sample} \times \text{Normality of FAS} \times 8000}{\text{Vol. of sample}}$$

Therefore COD of the water sample =

= .....500.....ppm

Water Analysis : Chemical Content



Copyright © Amrita University 2009 - 2015

**VARIABLES**

Select test:  
COD

Titrant :  
Domestic waste water (2...

Normality : 0.5

Speed of titrant : 0.2

START

RESET

**RESULT**

Titrant : Ferrous ammonium sulphate

Titrant used : 2.5 ml

Indicator : Ferroin

### COD of Tap-water

Volume of FAS used=  $(V_2 - V_1) = \dots 0.1 \dots \text{mL}$ .

Normality of FAS =  $\dots 0.5 \dots \text{N}$ .


Volume of the water sample =  $\dots 20 \dots \text{mL}$ .

$$\frac{\text{Vol. FAS for sample} \times \text{Normality of FAS} \times 8000}{\text{Vol. of sample}}$$

Therefore COD of the water sample =

=  $\dots 20 \dots \text{ppm}$

Water Analysis : Chemical Content



Copyright © Amrita University 2009 - 2015

**VARIABLES**

Select test:  
COD

Titrant:  
Tap water (20ml)

Normality : 0.5

Speed of titrant : 0.2

START  
RESET

**RESULT**

Titrant : Ferrous ammonium sulphate  
Titrant used : 4.9 ml  
Indicator : Ferroin

### Result:

COD of well-water sample =  $\dots 40 \dots \text{ppm}$ .

COD of Domestic-water sample =  $\dots 500 \dots \text{ppm}$ .

COD of tap-water sample =  $\dots 20 \dots \text{ppm}$ .